

セルカルチャーインサート適合表

-cell culture insert compatibility table-

Manufacturer	Part No.			Membrane				Insert Type
	24-Well Size	12-Well Size	6-Well Size	Material	Pore Size	Pore Density	Optics	
Corning Transwell	3413	3401	3412	PC	0.4 μm	1.0 × 10 ⁸ / cm ²	translucent	hanging
	3470	3460	3450	PET	0.4 μm	4.0 × 10 ⁶ / cm ²	transparent	hanging
	3415	3402	3414	PC	3.0 μm	2.0 × 10 ⁶ / cm ²	translucent	hanging
	3472	3462	3452	PET	3.0 μm	2.0 × 10 ⁶ / cm ²	transparent	hanging
Corning Falcon	353495	353494	353493	PET	0.4 μm	1.6 × 10 ⁸ / cm ²	translucent	hanging
	353095	353180	353090	PET	0.4 μm	1.6 × 10 ⁶ / cm ²	transparent	hanging
	353104	353103	353102	PET	1.0 μm	1.6 × 10 ⁶ / cm ²	transparent	hanging
	353492	353292	353092	PET	3.0 μm	2.0 × 10 ⁶ / cm ²	translucent	hanging
	353096	353181	353091	PET	3.0 μm	8.0 × 10 ⁵ / cm ²	transparent	hanging
Greiner Bio-One ThinCert	662 640	665 640	657 640	PET	0.4 μm	1.0 × 10 ⁸ / cm ²	translucent	hanging
	662 641	665 641	657 641	PET	0.4 μm	2.0 × 10 ⁶ / cm ²	transparent	hanging
	662 610	665 610	657 610	PET	1.0 μm	2.0 × 10 ⁶ / cm ²	transparent	hanging
	662 631	665 631	657 631	PET	3.0 μm	2.0 × 10 ⁶ / cm ²	translucent	hanging
	662 630	665 630	657 630	PET	3.0 μm	6.0 × 10 ⁵ / cm ²	transparent	hanging
Merck Millipore Millicell	MCHT 24H 48	MCHT 12H 48	MCHT 06H 48	PET	0.4 μm	1.0 × 10 ⁸ / cm ²	translucent	hanging
	MCRP 24H 48	MCRP 12H 48	MCRP 06H 48	PET	1.0 μm	2.0 × 10 ⁶ / cm ²	transparent	hanging
	MCSP 24H 48	MCSP 12H 48	MCSP 06H 48	PET	3.0 μm	2.0 × 10 ⁶ / cm ²	translucent	hanging
Brand	7827 00	7827 20	7827 40	PS	0.4 μm	1.0 × 10 ⁸ / cm ²	translucent	standing
Nunc	140620	140652	140660	PC	0.4 μm	8.5 × 10 ⁷ / cm ²	translucent	standing
	140627	140654	140663	PC	3.0 μm	1.7 × 10 ⁶ / cm ²	translucent	standing
Sarstedt	83.3932.040	83.3931.040	83.3930.040	PET	0.4 μm	1.0 × 10 ⁸ / cm ²	translucent	hanging
	83.3932.041	83.3931.041	83.3930.041	PET	0.4 μm	2.0 × 10 ⁶ / cm ²	transparent	hanging
	83.3932.101	83.3931.101	83.3930.101	PET	1.0 μm	2.0 × 10 ⁶ / cm ²	transparent	hanging
	83.3932.300	83.3931.300	83.3930.300	PET	3.0 μm	2.0 × 10 ⁶ / cm ²	translucent	hanging
SABEU cellQART	932 04 02	931 04 02	930 04 02	PET	0.4 μm	1.0 × 10 ⁸ / cm ²	transparent	hanging
	932 10 02	931 10 02	930 10 02	PET	1.0 μm	2.0 × 10 ⁶ / cm ²	transparent	hanging
	932 30 02	931 30 02	930 30 02	PET	3.0 μm	2.0 × 10 ⁶ / cm ²	transparent	hanging
	932 04 12	931 04 12	930 04 12	PET	0.4 μm	2.0 × 10 ⁶ / cm ²	transparent	hanging
	932 10 12	931 10 12	930 10 12	PET	1.0 μm	2.0 × 10 ⁶ / cm ²	transparent	hanging
	932 30 12	931 30 12	930 30 12	PET	3.0 μm	6.0 × 10 ⁵ / cm ²	transparent	hanging
Kanto chemical ad-MED Vitrigel 2	08364-96	08364-96	-	Collagen	-	-	translucent	hanging